

# NAP-Plasmid DNA

## 1. Why was there a decrease in the plasmid DNA recovery quantity?

Please confirm and follow the manual's instructions accordingly. There would not be a sufficient amount of plasmid if the lysis step was not performed thoroughly.

Elution buffer < pH 7.0, it was unable to confirm the eluted plasmid.

If the cultured bacterial species exhibited a low copy number of plasmid, the bacterial culture volume needed to be increased. Plasmid lost in E.coli, please re-culture with the freshly prepared bacteria. Elution buffer, pre-heat at 70°C.

## 2. Why was the RNA contamination present?

Please make sure RNase A was added into the Buffer M1 and stored at 4°C.

## 3. Why was the Inhibition of downstream enzymatic reactions present?

It indicated the residual presence of alcohol. After adding the wash Buffer2, re-centrifuge for the Centrifuge Step or place the column in the oven at 60°C for 10 minutes. For the RNA or DNA interference, remove and clean it well.

## 4. Why were there smeared plasmid bands on the agarose gel?

It indicated plasmid DNA degradation. Please shake it carefully when performing it. Keep plasmid preparations on ice.

## 5. Why was there poor plasmid quality?

Do not overgrow bacterial cultures. Do not incubate more than 5 minutes after adding the Buffer M1.