

# Ammonium persulfate , APS

Cat No. MB527-0010

Size: 10 g

## Description

The Ammonium persulfate (APS) is as a catalyst for acrylamide gel polymerization.

## Product Information

Synonym: AP, APS, Ammonium peroxodisulfate, Ammonium peroxydisulfate, PER

CAS number: 7727-54-0

EC index number: 016-060-00-6

EC number: 231-786-5

Grade: ACS, Reag. Ph Eur

Hill Formula:  $H_8N_2O_8S_2$

Chemical formula:  $(NH_4)_2S_2O_8$

Molar Mass: 228.19 g/mol

HS Code: 2833 40 00

## Application

Catalyst for acrylamide gel polymerization.

## Related Products

Cat. No.	Description	Size
MB526-0010	N,N,N',N'-Tetramethyl ethylenediamine (TEMED)	10 ml
NJ001-0010	Nimble juice Speedy Protein gel stain100X	10 ml
PM001-0500	BLUltra Prestained Protein Ladder	500 ul