



Media Formulation

0.05% Trypsin-EDTA, 1X Cat. No.: CC509-0100

Components	Concentration (mg/L)
Inorganic Salts	
Potassium Chloride (KCl)	400.0
Potassium Phosphate monobasic (KH ₂ PO ₄)	60.0
Sodium Bicarbonate (NaHCO ₃)	350.0
Sodium Chloride (NaCl)	8000.0
Sodium Phosphate dibasic (Na ₂ HPO ₄ ·7H ₂ O)	90.0
Other Components	
D-Glucose (Dextrose)	1000.0
EDTA·4Na·2H ₂ O	200.0
Phenol Red	10.0
Trypsin	500.0