

## 100 mM dNTP set

Part No./Catalog No.: DN040-4000

Lot Number: DN20340224

Minimum shelf life: Apr-2024

Size: 4x1ml

### Product Qualification

Specifications: The product must show the DNA amplification by PCR reaction. The bands when separated by 1.0% agarose electrophoresis with 1X TBE Buffer

Results: Meet specifications.

---

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



The DNA was amplified from p35S-GFP plasmid.

Gel photo taken under iBright FL1500 imaging system.

*Hsien yi Lin*

---

Hsien-Yi, Lin

Inspector, Quality

Issued on 2023, Apr