

100 mM dNTP set

Part No./Catalog No.: DN040-400M

Lot Number: DN21020252

Minimum shelf life: Dec-2022

Size: 4x100 ml

Product Qualification

Specifications: The product must show the DNA amplification by PCR reaction. The bands when separated by 1.0% agarose electrophoresis with 1X TBE Buffer

Results: Meet specifications.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



The DNA was amplified from p35S-GFP plasmid.

Gel was explored under the Blue light.

Hsien yi Lin

Hsien-Yi, Lin

Inspector, Quality

Issued on 2020, Dec