

## 100 mM dNTP set, PCR Grade

Part No./Catalog No.: DN046-1000

Lot Number: DN20350139

Minimum shelf life: May-2024

Size: 4x250ul

### Product Qualification

Specifications: The product must show the DNA amplification by PCR reaction. The bands when separated by 1.0% agarose electrophoresis with 1X TBE Buffer

Results: Meet specifications.

---

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



*Hsien yi Lin*

Hsien-Yi, Lin

Inspector, Quality

Issued on 2023, May

The DNA was amplified from p35S-GFP plasmid.

Gel photo taken under iBright FL1500 imaging system.