

Genomic DNA Isolation Kit (Tissue)

Catalog Number: NA026-0100

Lot Number: 70101160

Size: 100 reactions

Quality Control

Specifications:

Content	TL	TP	W1	W2	BE	TC Column	Collection Tube
Lot.No	70101161	70101162	70101163	70101164	70101165	70101166	70101167
pH	8.5	5.5	8.5	8.5	8.5	----	----
accredited	√	√	√	√	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Pork	30 mg	21.3	√	√
2	Pork	30 mg	21.6	√	√
3	Pork	30 mg	22.1	√	√
4	Pork	30 mg	22.1	√	√
5	Pork	30 mg	21.8	√	√
6	Pork	30 mg	22.1	√	√

Result: Meets specification.

Storage condition: Genomic DNA Isolation Kit (Tissue) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.