

Genomic DNA Isolation Kit (Tissue)

Catalog Number: NA026-0100

Lot Number: 70109150

Size: 100 reactions

Quality Control

Specifications:

Content	TL	TP	W1	W2	BE	TC Column	Collection Tube
Lot.No	70109151	70109152	7010913	70109154	70109155	70109156	70109157
pH	8.5	5.5	8.5	8.5	8.5	----	----
accredited	√	√	√	√	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Pork	30 mg	21.3	√	√
2	Pork	30 mg	21.6	√	√
3	Pork	30 mg	22.1	√	√
4	Pork	30 mg	22.1	√	√
5	Pork	30 mg	21.8	√	√
6	Pork	30 mg	22.1	√	√

Result: Meets specification.

Storage condition: Genomic DNA Isolation Kit (Tissue) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.