

50mL Serological Pipet



Cat. No. PC550-0100

Size: 100 Pipets/Case

Description: 50mL in 5/10mL, Individually wrapped, paper / plastic

Sterile: Yes


Purpose: Use for liquid transfer

Materials

Pipet: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Color: Clear



Unit: mm

	Dimensions	Tolerance
A	∅18.4	±0.10
B	∅ 8.0	±0.25
C	∅ 5.8	±0.25
D	346.6	±2.0

Features

- Temperature: -20°C to +50°C
- Sterilized by gamma irradiation SAL10⁶(ISO11137)
- Latex Statement: this product is latex free
- DNase/RNase free, non-pyrogenic, non-cytotoxic
- Meet the requirements of BSE/TSE
- Use dispenser box for easy opening
- Available in Individually packed in peel-to-open wrap
- Autoclavability: No
- Negative graduations up to -10 ml (total volume = 60 ml)
- Graduation (mL): in 5/10mL
- Graduations are calibrated for accurate dispensing to within ±2%
- Color coded: purple
- Shelf Life: 3 years after month of production
- Every product is printed with a batch number for traceability
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system