

Certificate of Analysis



Novel Juice

Supplied in 6X Loading Buffer

Part No./Catalog No. SL001-0050

Lot Number: SL20340240

Minimum shelf life: Apr-2025

size: 50µl

1 ul per well for gel applications


Product Qualification

Specifications: In 1% agarose gels, xylene cyanol FF co-migrates with ~4000 bp DNA, while bromophenol blue co-migrates with ~300 bp DNA, and Orange G co-migrates with smallest fragments, ~50 bp.

Results: Meet specifications.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.



Hsien yi Lin

Hsien-Yi, Lin

Inspector, Quality

Issued on 2023, Apr

1ul SL001-0050, with 5ul 100bp DNA Ladder.

Gel photo taken under iBright FL1500 imaging system.