

Genomic DNA Isolation Kit (Tissue)

Catalog Number: SN026-0100

Size: 100 reactions

Catalog Number: SN026-0004

Size: 4 reactions

Lot Number: 90109190

Quality Control

Specifications:

We show the product pH and recovery ratio in below:

Content	TL	TP	W1	W2	BE	TC Column	Collection Tube
Lot No	90109191	90109192	90109193	90109194	90109195	90109196	90109197
pH	8.5	5.5	8.5	8.5	8.5	----	----
accredited	√	√	√	√	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Pork	30 mg	21.0	√	√
2	Pork	30 mg	22.1	√	√
3	Pork	30 mg	22.3	√	√
4	Pork	30 mg	22.1	√	√
5	Pork	30 mg	21.5	√	√
6	Pork	30 mg	22.4	√	√

Result: Meets specification.

Storage condition: Genomic DNA Isolation Kit (Tissue) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.