

Certificate of Analysis



DMEM, High Glucose

Catalog Number: CC103-0500

Lot Number: GL20001

Size: 500 mL

Minimum Shelf life: Feb-2021

Storage: 2 ~ 8 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.04 EU/mL
Osmolality	Check & Record	344
pH	Check & Record	7.3
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.