

# Certificate of Analysis



## DMEM, High Glucose

Catalog Number: CC107-0500

Lot Number: GL20099

Size: 500 mL

Minimum Shelf Life: Jul-2021

Storage: 2 ~ 8 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.04 EU/mL
Osmolality	Check & Record	312
pH	Check & Record	7.1
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

A handwritten signature in black ink that reads 'Yi Hsing, Yao'.

Yi-Hsing, Yao  
Inspector, Quality  
Issued on 2020/07/30