

Certificate of Analysis



DMEM Powder, High Glucose

Catalog Number: CC126-1010

Lot Number: GL20128

Size: 10 Liter

Minimum Shelf Life: Sep-2023

Storage: 2 ~ 8 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.004 EU/mL
Osmolality	Check & Record	266
pH	Check & Record	6.2
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Yi-Hsing, Yao
Inspector, Quality
Issued on 2020/10/07