

Certificate of Analysis



TTQ12140



PBS (10X), pH 7.4

Catalog Number: CC708-0500

Lot Number: GL20095

Size: 500 mL

Minimum Shelf Life: Jul-2022

Storage: 2 ~30 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.08 EU/mL
Osmolality	Check & Record	294
pH	Check & Record	7.4
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Yi-Hsing, Yao

Inspector, Quality

Issued on 2020/07/30