

Certificate of Analysis



TTQ12140



IT-101578

10X PBST

Catalog Number: CC710-0500

Lot Number: GL20096

Size: 500 mL

Minimum Shelf Life: Jul-2022

Storage: 2 ~30 °C

Test	Specifications	Results
Endotoxic Assay	Check & Record	<0.04 EU/mL
Osmolality	Check & Record	455
pH	Check & Record	7.5
Sterility Testing	Negative	Negative

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Yi-Hsing, Yao

Inspector, Quality

Issued on 2020/07/30