

Certificate of Analysis

XLarge DNA Ladder RTU

Part No./Catalog No. DM013-R025

Lot Number: DM11520032

Size: 25ul


Product Qualification

Specifications: The product must show 14 dsDNA bands when separated by electrophoresis on a 0.7% agarose gel with 1x TBE buffer.

Results: Meet specifications.

Specification: Any contaminant band must be less than 5 % in intensity to the nearest major band shown by EtBr staining.

Result: Meets specification.



All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

www.genedirex.com
info@genedirex.com

GeneDireX[®]