

Certificate of Analysis

XLarge DNA Ladder RTU

Part No./Catalog No. DM013-R500

Lot Number: 11DM504102

Size: 50µg / 500µl

Product Qualification

Specifications: The product must show 14 dsDNA bands when separated by electrophoresis on a 0.7% agarose gel with 1x TBE buffer.

Results: Meet specifications.

Specification: Any contaminant band must be less than 5 % in intensity to the nearest major band shown by EtBr staining.

Result: Meets specification.



All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.