

# Certificate of Analysis

## OneMARK B DNA Ladder

Part No./Catalog No. DM110-0100

Lot Number: 13DM032701

Size: 50µg/ 600µl

6 ul per well for 100 gel applications

### Product Qualification

Specifications: The product must show 12 dsDNA bands when separated by electrophoresis on a 1% agarose gel with 1x TBE buffer.

Results: Meet specifications.


Specification: Any contaminant band must be less than 5 % in intensity to the nearest major band shown by Novel Juice staining.

Result: Meets specification.

---

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.



6ul , 1% agarose gel with 1x TBE buffer