

Certificate of Analysis



GScript RTase

Part No./Catalog No. MB303-0050

Lot Number: MB20361107

Minimum shelf life: Jun-2024

Size: 50 Reactions

Product Qualification

Specifications:

GScript RTase is tested functionally using 293cell RNA as the template for amplification of a 1176bp target (30 cycles).

This product was tested for amplification of control gene from 293cell cDNA.

Results: Meet specifications.



Hsien yi Lin

Hsien-Yi, Lin

Inspector, Quality

Issued on 2023, Jun

Human cell cDNA amplification product

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.