

Certificate of Analysis



Random Hexamer Primer

Part No./Catalog No. MB702-0100

Lot Number: MB20040221

Concentration: 50uM


Product Qualification

Specifications:

RT-PCR System is tested functionally using 293cell RNA as the template for amplification of a 1176bp target (25 cycles).

This product was tested for amplification of 1176 bp gene from 293cell cDNA.

Results: Meet specifications.



Hsien-Yi Lin

Hsien-Yi, Lin
Inspector, Quality
Issued on 2020/04

Lane M: 1kb DNA Ladder (GeneDirex, DM010-R500)

Lane 1: 293cell cDNA amplification product

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic **USES**.