

Certificate of Analysis



EVOgel

Part Number/Catalog Number: MB803-0100

Lot Number: 90104080

Size: 100 ml

Expiry Date: Jun-2021

Product Qualification

Appearance: Clear, colorless liquid.

Mol. Formula: Not Applicable

Mol. Weight: Not Applicable

Description:

Assay	Min	Max	Units	Disposition
Conductivity (1:20, water)	430	550	US	PASS
Electrophoresis	PASS		P/F	PASS
pH @25°C	6.8	7.2		PASS
Polymerization Time		30	minutes	PASS
Protease		NONE	P/F	PASS

Storage: Store at room temperature. Protect from moisture.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.