

# Certificate of Analysis

## PG Reagent

### (Plant Genomic DNA Isolation Kit)

Part No./Catalog No.: NA002-0004

Lot Number: 12NA040200

Size: 4 reactions

#### Quality Control

Content	PG
pH	5.2
accredited	√

Result	Test Sample	Test Volumn	Yield (µg)	Agarose check	accredited
1	Orchid leaf	100 mg	14	√	√
2	Orchid leaf	100 mg	15	√	√
3	Orchid leaf	100 mg	16	√	√
4	Orchid leaf	100 mg	15	√	√
5	Orchid leaf	100 mg	16	√	√
6	Orchid leaf	100 mg	16	√	√

#### Storage conditions

PG Reagent should be stored dry at room temperature (15–25°C).

---

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses