

## MBead Bacteria Genomic DNA Kit

Catalog Number: NA010-0100

Lot Number: 13NA011600

Size: 100 reactions

## Product Qualification

Specifications:

We show the product pH and recovery ratio in below:

Content	Lysis Buffer	Wash Buffer	Release Buffer	Magnetic Bead
pH	4.5	5.0	8.5	---
accredited	√	√	√	√

Result	Test Sample	Test Volumn	Yield (µg)	Agarose check	accredited
1	DH5α	300µ l	31.9	√	√
2	DH5α	300µ l	32.9	√	√
3	DH5α	300µ l	33.8	√	√
4	DH5α	300µ l	30.4	√	√
5	DH5α	300µ l	33.2	√	√
6	DH5α	300µ l	31.8	√	√

Result: Meets specification.

**Storage condition:** MBead Bacteria Genomic DNA Kit should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses