

Certificate of Analysis

MBead Plant Genomic DNA Kit

Catalog Number: NA012-0004 Size: 4 reactions

Lot Number: 301122400

Product Qualification

Specifications:

We show the product pH and recovery ratio in below:

Content	Lysis Buffer	Grind Buffer	Wash Buffer	Release Buffer	Magnetic Bead
Lot.No.	301122402	301122403	301122404	301122405	301122401
pH	4.5	8.5	5.0	8.5	---
accredited	√	√	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Orchid leaf	100 mg	15.7	√	√
2	Orchid leaf	100 mg	15.5	√	√
3	Orchid leaf	100 mg	16.8	√	√
4	Orchid leaf	100 mg	14.9	√	√
5	Orchid leaf	100 mg	15.1	√	√
6	Orchid leaf	100 mg	15.4	√	√

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses