

MBead Blood/Cell Genomic DNA Kit

Catalog Number: NA013-0100

Lot Number: 12NA122100

Size: 100 reactions

Quality Control

Content	Magnetic Bead	Lysis Buffer	Wash Buffer	Release Buffer
pH	4.5	5.0	8.5	---
accredited	√	√	√	√

Result	Test DNA Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Whole Blood	300µ l	12.9	√	√
2	Whole Blood	300µ l	13.8	√	√
3	Whole Blood	300µ l	11.8	√	√
4	Whole Blood	300µ l	12.1	√	√
5	Whole Blood	300µ l	11.2	√	√
6	Whole Blood	300µ l	14.1	√	√

Storage conditions

MBead Blood/Cell Genomic DNA Kit should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses