

Certificate of Analysis



Total RNA Isolation Kit (Plant)

Catalog Number: NA020-0004

Size: 4 reactions

Lot Number: 70105110

Product Qualification

Specifications:

We show the product pH and recovery ratio in below:

Content	RP	W1	W2	RE	RP Column	Collection Tube
Lot.No.	70105111	70105112	70105113	70105114	70105115	----
pH	5.6	8.5	8.5	8.5	----	----
accredited	√	√	√	√	√	√

Result	Test RNA Sample	Test	Yield (µg)	Agarose check	accredited
1	Orchid leaf	100 mg	5.5	√	√
2	Orchid leaf	100 mg	5.6	√	√
3	Orchid leaf	100 mg	5.6	√	√

Result: Meets specification.

Storage condition: Total RNA Isolation Kit (Plant) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses