

## Total RNA Isolation Kit (Plant)

Catalog Number: NA020-0100

Lot Number: 70101160

Size: 100 reactions

## Product Qualification

Specifications:

We show the product pH and recovery ratio in below:

Content	RP	W1	W2	RE	RP Column	Collection Tube
Lot.No.	70101161	70101162	70101163	70101164	70101165	70101166
pH	5.6	8.5	8.5	8.5	----	----
accredited	√	√	√	√	√	√

Result	Test RNA Sample	Test	Yield (µg)	Agarose check	accredited
1	Orchid leaf	100 mg	5.5	√	√
2	Orchid leaf	100 mg	5.6	√	√
3	Orchid leaf	100 mg	5.6	√	√

Result: Meets specification.

**Storage condition:** Total RNA Isolation Kit (Plant) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses