

15mL Conical Centrifuge Tube

Cat. No. PC116-0500

Size: 50/Pack, 500/Case



Description

Temperature Range: Approx.-80 to+121
Shelf Life: 3 years after date of production
Sterile: Yes

Product Range

Volume(ML): 15 ML
Color: transparency
Packaging: Bagged
Max.RCF(xg): 12,000

Features

- USP VI Materials
- Meet the requirements of analytical level
- High strength non-toxic materials
- Printed clear graduations with marking area
- PP for the tube and HDPE for the cap
- One-hand operation design
- Sterilized by E-Beam
- DNase, RNase and pyrogen free

Technical Drawing

