

100mm Cell Culture Dish



Cat. No. PC205-0300

Size: 300 Pipets/Case

Spec		Cell Growth Area (cm ²)	Recommended Medium Volume (mL)
Type	Height (mm)		
100	20	57.6	11-16.5

Description

Temperature Range:

Approx. -86 to +64°C

Shelf Life:

3 Years after date of production (ensure package is in good)

Sterile:

Yes

Features

- High clarity polystyrene
- Flat transparent surface for distortion-free observation
- Vacuum plasma TC treatment, excellent cell adherence
- Stackable for easy storage and handling
- Sterilized by E-Beam
- Non-pyrogenic

