

# 5mL Serological Pipet

Cat. No. PC505-0200

Size: 200 Pipets/Case



**Description:** 5mL in 1/10mL, Individually wrapped, paper / plastic

**Sterile:** Yes


**Purpose:** Use for liquid transfer

## Materials

Pipet: GPPS (General polystyrene), Compliance with USP Class VI

Filter: PO (Polyolefin)

Color: Clear



Unit: mm

	Dimensions	Tolerance
A	Ø8.2	±0.1
B	Ø2.5	±0.35
C	341.0	±2.0

## Features

- Temperature: -20°C to +50°C
- Sterilized by gamma irradiation SAL10<sup>6</sup>(ISO11137)
- Latex Statement: this product is latex free
- DNase/RNase free, non-pyrogenic, non-cytotoxic
- Meet the requirements of BSE/TSE
- Use dispenser box for easy opening
- Available in Individually packed in peel-to-open wrap
- Autoclavability: No
- Negative graduations up to -3 ml (total volume = 8 ml)
- Graduation (mL): in 1/10mL
- Graduations are calibrated for accurate dispensing to within ±2%
- Color coded: blue
- Shelf Life: 3 years after month of production
- Every product is printed with a batch number for traceability
- Manufactured in a class 100,000 room environment
- Manufactured under ISO13485 and ISO9001 quality management system