

BlueRAY Prestained Protein Ladder

Part No./Catalog No.: PM006-0500

Lot Number: 70251009

Size: 500 ul

Product Qualification

Specifications: The product must show 10 major protein bands when separated by SDS-PAGE electrophoresis on a 15% Tris-Glycine Gel with Tris-Glycine buffer, 12% Bis-Tris Gel with MOPS buffer, and 4-12% Bis-Tris Gel with MES buffer.

Results: Meet specifications.

Specification: Any pre-stained contaminant band must be less than 10% in intensity to the nearest major band.

Result: Meets specification.

Specification: Pre-stained protein band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic USES.



5ul, 4-20% Tris- Glycine Gel with Tris-Glycine buffer