

BlueRAY Prestained Protein Ladder

Catalog Number: PM006-500M

Size: 500 ml

Lot Number: 91068021

Product Qualification

Specifications: The product must show 10 major protein bands when separated by SDS-PAGE electrophoresis on a 15% Tris-Glycine Gel with Tris-Glycine buffer, 12% Bis-Tris Gel with MOPS buffer, and 4-12% Bis-Tris Gel with MES buffer.


Results: Meet specifications.

Specification: Any pre-stained contaminant band must be less than 10% in intensity to the nearest major band.

Result: Meets specification.

Specification: Pre-stained protein band from lot being tested must migrate within 1/2 bandwidth from prior lot.

Result: Meets specification.



5 μ l, EVOGel with Tris-Glycine-SDS running buffer

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.