



Certificate of Analysis

Genomic DNA Isolation Kit

(Blood/Cultured Cell/Fungus)

Catalog Number: SN023-0004

Size: 4 reactions

Lot Number: 70112190

Quality Control

Specifications:

Content	CR	CC	CB	W1	W2	BE	CC Column	Collection Tube
Lot No.	70112191	70112192	70112193	70112194	70112195	70112196	70112197	----
pH	8.1	8.5	4.5	8.5	8.5	8.5	----	----
accredited	√	√	√	√	√	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Whole Blood	200µl	5.2	√	√
2	Whole Blood	200µl	5.2	√	√
3	Whole Blood	200µl	5.4	√	√
4	Whole Blood	200µl	5.5	√	√
5	Whole Blood	200µl	5.7	√	√
6	Whole Blood	200µl	5.6	√	√

Result: Meets specification.

Storage condition: Genomic DNA Isolation Kit (Blood/Cultured Cell/Fungus) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.