

Genomic DNA Isolation Kit (Blood/Cultured Cell/Fungus)

Catalog Number: SN023-0100 Size: 100 reactions

Catalog Number: SN023-0004 Size: 4 reactions

Lot Number: 00102040

Minimum Shelf Life: Feb-2022

Quality Control

Specifications:

Content	CR	CC	CB	W1	W2
Lot No.	00102041	00102042	00102043	00102044	00102045
pH	8.1	8.5	4.5	8.5	8.5
accredited	√	√	√	√	√

Content	BE	CC Column	Collection Tube
Lot No.	00102046	00102047	00102048
pH	8.5	----	----
accredited	√	√	√

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	accredited
1	Whole Blood	200µl	5.0	√	√
2	Whole Blood	200µl	5.1	√	√
3	Whole Blood	200µl	5.5	√	√
4	Whole Blood	200µl	5.4	√	√
5	Whole Blood	200µl	5.5	√	√
6	Whole Blood	200µl	5.0	√	√

Result: Meets specification.

Storage condition: Genomic DNA Isolation Kit (Blood/Cultured Cell/Fungus) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Certificate of Analysis



Caution: Not intended for human or animal diagnostic or therapeutic uses.

Yi Hsing, Yao

Yi-Hsing, Yao

Inspector, Quality

Issued on 2020/ Feb/04