

## Genomic DNA Isolation Kit (Blood/Cultured Cell/Fungus)

Catalog Number: SN023-0100      Size: 100 reactions  
 Catalog Number: SN023-0004      Size: 4 reactions  
 Lot Number: 00112180

### Quality Control

Specifications:

Content	CR	CC	CB	W1	W2
Lot No.	00112181	00112182	00112183	00112184	00112185
pH	8.1	8.5	4.5	8.5	8.5
accredited	√	√	√	√	√

Content	BE	CC Column	Collection Tube
Lot No.	00112186	00112187	00112188
pH	8.5	----	----
accredited	√	√	√

*Yi Hsing, Yao*

Yi-Hsing, Yao  
 Inspector, Quality  
 Issued on 2020/ 12/18

Result	Test Sample	Test Volume	Yield (µg)	Agarose check	Accredited
1	Whole Blood	200µl	5.7	√	√
2	Whole Blood	200µl	5.8	√	√
3	Whole Blood	200µl	5.3	√	√
4	Whole Blood	200µl	5.9	√	√
5	Whole Blood	200µl	5.6	√	√
6	Whole Blood	200µl	5.9	√	√

Result: Meets specification.

**Storage condition:** Genomic DNA Isolation Kit (Blood/Cultured Cell/Fungus) should be stored dry at room temperature (15–25°C).

All products are for research use only.

Caution: Not intended for human or animal diagnostic or therapeutic uses.