

225cm² Cell Culture Flask



Cat. No.	Cell Growth Area(cm ²)	Cap Style	Volume (mL)	Recommended Medium Volume(mL)	/Pack	/Case
PC279-0025	225	Vent Cap	950	45-67.5	5	25

Simply cell culture flasks are an ideal solution for enhanced cell attachment and growth of a variety of primary cells and transformed cells in serum free or serum-containing cultures.

Product Range

Temperature Range: Approx. -86 to +64°C

Shelf Life: 3 Years after date of production (ensure package is in good)

Sterile: Yes

Features

- High clarity medical grade polystyrene material
- Sterilized by E Beam SAL=10⁻⁶
- Non-pyrogenic, DNase/ RNase free, free, Non-Cytotoxicity
- Packaged in sterile, zip sealable bags
- Two types of cap style: vent caps
- Caps with 0.22um hydrophobic filters for gas exchange without contamination
- Fostered writing and clear graduation
- Clear lot number for batch traceability
- Stackable

