

60mm Cell Culture Dish



Cat. No. PC202-0500

Size: 500 Dishes/Case

Spec		Cell Growth Area (cm ²)	Recommended Medium Volume (mL)
Type	Height (mm)		
60	15	22.9	4.2-6.3

Description

Temperature Range:

Approx. -86 to +64°C

Shelf Life:

3 Years after date of production (ensure package is in good)

Sterile:

Yes

Features

- High clarity polystyrene
- Flat transparent surface for distortion-free observation
- Vacuum plasma TC treatment, excellent cell adherence
- Stackable for easy storage and handling
- Sterilized by E-Beam
- Non-pyrogenic

