

60mm Cell Culture Dish (Gripping Ring)



Cat. No. PC203-0600

Size: 600 Dishes/Case

Spec		Cell Growth Area (cm ²)	Recommended Medium Volume (mL)
Type	Height (mm)		
60 (With Gripping Ring)	14.78	21.2	4-7

Description

Temperature Range:

Approx. -20°C~+50°C

Shelf Life:

3 years after month of production

Sterile:

Yes

Features

- Every inner bag is printed with a batch number for traceability
- Sterilized by gamma irradiation SAL10⁻⁶ (ISO11137)
- Non-pyrogenic
- Flat bottom uniform wall thickness ensures distortion-free bottom
- Lids with several little chimbs to shape vents are available for very effective gas exchange
- Bulk packed in zip-Lock bag
- Autoclavability: No

- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485 and ISO9001 quality management system

Dimensions (Unit: mm)

